



## Print Recipe

Baby Boomer Gel Rose

CPNP 281 8071

Typ	Deklaration	Cas-Nr.	FDA
INCI	DI-HEMA TRIMETHYLHEXYL DICARBAMATE	72869-86-4	A2
INCI	URETHANE ACRYLATES/METHACRYLATES	n.b.	B
INCI	HYDROXYPROPYL METHACRYLATE	27813-02-1	D
INCI	URETHANE ACRYLATES	n.b.	D
INCI	GLASS	65997-17-3	E
INCI	SILICA DIMETHYL SILYLATE	68611-44-9	E
INCI	TRIPROPYLENE GLYCOL DIACRYLATE	42978-66-5	E
INCI	METHYL BENZOYLFORMATE	15206-55-0	F
INCI	BHT	128-37-0	G
INCI	STEARIC ACID	57-11-4	G
INCI	ALUMINA	1344-28-1	G
INCI	P-HYDROXYANISOLE	150-76-5	G
INCI	CI 77891 (Nano)	13463-67-7	F
INCI	CI 77891	13463-67-7	G
INCI	CI 15985:1	2783-94-0	G
INCI	CI 45410:2	18472-87-2	G
INCI	CI 15850:2	5858-81-1	G
INCI	CI 77266	1333-86-4	G

FDA: A1 = 75–100 %; A2 = 50-75 %; B = 25-50%, C = 10-25%; D = 5-10%; E = 1-5%; F = 0,1-1%; G = 0-0,1%, H = Traces

### Product description:

Medium-viscosity gel in a pink shade that, when combined with Baby Boomer Gel 'white', makes it quick and easy to achieve the perfect transition in a Baby Boomer design. Excellent adhesion properties and no yellowing.

### Application instructions:

Apply the gel to the prepared natural fingernail or tips in one or two thin layers. The curing time per layer is 60-90 seconds under UV/LED light.

## Baby Boomer Gel Rose

### Labelling Elements:



Danger

H317 – May cause allergic skin reactions

H319 – Causes serious eye irritation

H411 – Toxic to aquatic life with long-lasting effects

EUH204 – Contains isocyanates. May cause allergic reactions.

P280 – Wear protective gloves made of butyl rubber in accordance with EN 374 and eye protection in accordance with EN 166

P302 + P352 – IF ON SKIN: Wash with plenty of soap and water

P333 + P313 – If skin irritation or rash occurs: Seek medical advice/attention.

P305 + P352 + P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337 + P313 – If eye irritation persists: Get medical advice/attention

P362 + P364 – Take off contaminated clothing and wash before reuse.

P273 – Avoid release to the environment.

P501 – Dispose of contents/container to an approved waste disposal company.

PAO – If stored correctly, can be kept for 12 months after opening (can symbol).

Unopened in its original packaging and stored correctly (cool, dry, protected from sunlight between 10 and 25°C), can be kept for several years.

UN number 3082